



*Feeling is Believing...
Experience the LIFT™ Magic...p.6*

*Check our website to see
our competitive
prices!*

House of Acupuncture

danish distributor of



www.houseofacupuncture.dk

House of Acupuncture

House of Acupuncture is the new serious supplier for the Danish TCM community.

To prove our serious way of business, we sell the well-proven ISO and CE certified products from HELIO, one of the leading manufacturer of quality needles and supplies for acupuncture.

HELIO sees it as their responsibility to prove you with the finest quality products and service. In 1997, they restructured the company and implemented the US FDA QSR and ISO quality systems. In that same year, they received both ISO and EN-46002 certificates as well as CE mark certification.

In 2006 they upgraded their quality system to ISO 13485:2003 with CMDCAS accreditation. To date, HELIO remains the only ISO certified company in the U.S. TCM supply industry. Their manufacturing sites where they produce the quality certified products are ISO 13485:2003 certified as well.

House of Acupuncture is committed to provide quality products, prompt delivery and reliable service. Together with HELIO and the Danish TCM community, we are always working to improve our product assortment and our service. As we say: "no one is better than their reputation" so always feel free to contact House of Acupuncture if you have any questions or ideas at all.

Last but not least, we are very proud to serve the Danish TCM community with the high quality products from HELIO.



*Now you can shop online.
Ordering supplies just got easier!*
www.houseofacupuncture.dk

Note: Prices subject to change without prior notice.

Contents

About the company
About The Company 2

Acupuncture Needles
Introducing "LIFT™" 6
Acupuncture Needles Overview 7
Innovative Needle Packaging Design 8
AcuGlide® 9
VINCO® 11
C&G™ 13

7-Star Needle, Press Tacks, and Pellets
Seven Star & Intradermal Needles 14
Press Tacks 15
Pellets 16

Moxibustion
Maxim™ Moxa 17
Moxa Sticks 18
Needle Moxa, and Loose Moxa 19
Stick-On Moxa Cones 20
Moxa Tools 21

Clinical Supplies
Needle Disposal Containers 22
Diagnostic Tools 22
Medical Supplies 22
Clinical Supplies 23

Treatment Devices
Cupping Devices 24

Therapeutic Equipment
Warranty Information 25
Stimulators and Point Detector 26
Electro-Acupuncture Devices 27
TENS 29
Massage Tables, Chairs and Accessories 32

Models and Charts
Models and Charts 35



For Our Supplier Helio, Quality Is A Way Of Life

Few companies invest in quality the way we do at Helio. In fact, our total commitment to quality products and service drives every decision we make. Helio's goal is for all our products to meet your needs and exceed your expectations. We are dedicated to maintaining our role as the undisputed leader in overall quality. At Helio we never rest in our search for advances in quality consistency and manufacturing excellence that hold the promise of further improving our products and service. Other companies might trumpet their claims of high quality products, but Helio backs up its quality assurance with significant investments in manufacturing and quality control facilities that are unsurpassed in the industry.



Our Supplier Helio Delivers Quality, Innovation, and Value!

Helio's electronic product testing lab performs burn-in and submicro-ampere analysis for our electronic products.

We contract with third party labs to conduct critical studies and tests including bio-compatibility tests, scanning micrograph analysis, x-ray dispersion analysis, infrared Fourier transform analysis, etc.

Helio also contracts with independent biological testing laboratories to validate our needle package designs for microbiological parameters, such as the biocompatibility of the needle bodies, microbial barrier effective period evaluation and pyrogen residual analysis.

Suzhou Helio Tech Co. located at Suzhou, China is our wholly-owned subsidiary - the HTC is fully-equipped with advance manufacturing and test equipment dedicated to the fabrication of quality acupuncture needles and product development.



Needle blister package under integrity evaluation




Needle inspection with digital microscopes



Simulated needle insertion smoothness test



The FIRARD-3 burn-in line up



We Try Hard to Earn Your Trust and Confidence.




At Helio, Quality is A Way of Life!



Acupuncture Needle Redefined!

LIFTed
Micro-Channel
Needle Surface



SEM at 300X



Typical
Smooth
Needle Surface



SEM at 300X

The Low Insertion Force Technology

You can feel it. Patient appreciates it.

WHAT IS LIFT™?

LIFT™ is the acronym of **Low Insertion Force Technology**. It is a Helio invention with patent pending. LIFT™ allows you to insert a non-silicone coated acupuncture needle with significantly less force compared to a typical acupuncture needle. It is beyond smooth.

WHAT IS THE BASIS OF LIFT™?

The LIFTed™ needle has a MicroChannel™ surface structure that significantly reduces the friction between the needle and tissue cells. Hence, it significantly reduces the forces required to insert the needle.

WHY USE LIFTED™ NEEDLE?

When you use a LIFTed™ acupuncture needle, you can immediately feel an effortless insertion. For your patient, he may not even feel being needled.

NOTE: House of Acupuncture offers LIFTed™ acupuncture needles in different brands, which are the best of its kind. Choose from AcuGlide®, AcuMaster®, VINCO®, and VIVA®.

Acupuncture Needles Overview

CE 0086

Caution:
Federal law restricts this device (acupuncture needles) to sale by or on the order of qualified practitioners of acupuncture as determined by the states.

Key Features:

- Surgical stainless steel needle body
- LIFTed™ needle body
- Bio-safe surface treatment
- Light weight aluminum alloy handle
- Copper handle
- SkinSaver™ guide tube
- FingerSaver™ packet
- Single needle with guide tube in blister pack
- Single needle without guide tube in blister pack
- ZiPack®
- 5-needle cluster pack
- EZY-Pack®
- Detox in Q-Pack™

Applications:

	 AcuGlide®	 vinco™ MicroClean®	 C&G™
Key Features:			
• Surgical stainless steel needle body	Yes	Yes	Yes
• LIFTed™ needle body	Yes	Yes	X
• Bio-safe surface treatment	Yes	X	X
• Light weight aluminum alloy handle	Yes	Yes	X
• Copper handle	X	X	Yes
• SkinSaver™ guide tube	Yes	Yes	X
• FingerSaver™ packet	Yes	Yes	X
• Single needle with guide tube in blister pack	Yes	Yes	Yes
• Single needle without guide tube in blister pack	Yes	Yes	Yes
• ZiPack®	Yes	X	X
• 5-needle cluster pack	X	Yes	Yes
• EZY-Pack®	Yes	Yes	X
• Detox in Q-Pack™	X	Yes	X
Applications:			
• Japanese-style acupuncture	• Japanese-style acupuncture	• General acupuncture	• General acupuncture
• When painless insertion is a must for extremely sensitive patients	• When painless insertion is a must for extremely sensitive patients	• Quick needle access	• When concerned about silicone coating
• When concerned about silicone coating	• When concerned about silicone coating	• Detoxification treatment	• When concerned about silicone coating
• When nickel allergic reaction is a concern	• When nickel allergic reaction is a concern	• When concerned about silicone coating	• When needle cost is a major concern
• Thin needle for children	• Thin needle for children		
• Thin needle for facial treatment	• Thin needle for facial treatment		
• When patients demand the best	• When patients demand the best		

Innovative Needle Package Design

For greater patient comfort and enhanced practice efficiency!

The following package designs were developed to enhance your patient's perception of your practice, save your valuable time, protect your fingers, and reduce clutter in your treatment room:

SkinSaver™ Guide Tube:



The traditional guide tube has sharp edges. The sharp edges scratch the patient's skin and cause pain, particularly when inserting the needle at an angle. This pain may create the perception in a patient that acupuncture is painful or that your needling technique is poor.

The SkinSaver™ guide tube has rounded edges, which eliminates the possibility of scratching the patient's skin. SkinSaver is available in AcuGlide®, AcuMaster® and VINCO™, except for the Cluster Pack.

FingerSaver™ Packet:



Most needle packages use stiff plastic blister which has sharp edges. While using stiff plastic may have its own advantages, the sharp edges present a hazard to the practitioner's fingers, particularly when trying to obtain a needle in a hurry.


The FingerSaver packet has rounded edges at the opening edge of the packet to eliminate the pricks and scratches to your fingers! FingerSaver is available in AcuGlide®, AcuMaster® and VINCO™.

ZiPack®

This package design is for needles with guide tubes, but with no tabs or glue spots for affixing the needle to the guide tube. Not having to clean the tabs off the floor, or struggle with sub-glued needles will save you time and money.

The ZiPack® is available in AcuGlide®, AcuMaster®.

Easy Step to Maximize the Advantage of ZiPack®




Innovatively Engineered Fine Needle with Extraordinary Features...

Key Features:

- LIFTed™ needle surface for maximum patient comfort and ease of insertion.
- Surface treated with parylene, no residual left behind.
- BioSafe™ – Class VI compatibility meeting US FDA G95-1 and ISO-10993 requirements as substantiated by third party GLP laboratory.
- Non-tarnishing aluminum pipe handle as light as plastic handles.
- Consistently sharp needle tip profile.
- SkinSaver™ painless guide tube for patient comfort.
- FingerSaver™ prick-free packet.
- Smallest guide tube for improved insertion precision.
- EZY-Accessible box design.
- Available in EZY-Pack® for busy practitioners.
- ZiPack™ – No tabs to mess with, no glue spot to fight with, just aim and insert!
- Size color-coded box and packet.

Applications:


- Japanese-style acupuncture
- When painless insertion is a must for extremely sensitive patients
- When nickel allergic reaction is a concern
- Thin needle for children
- When concerned about silicone coating
- Thin needle for facial treatment
- When patients demand the best



CE 0086



The Fine Needle with Extraordinary Features...



MT1: ZiPack®



- Individually packaged with SkinSaver™ guide tubes.
- No glue spots.
- No tabs.
- Aim and insert.
- 5 needles per flat with FingerSaver™ feature.
- 20 flats; 100 needles per box.




ML1: TL-Pack®



- Individually packaged without guide tube in FingerSaver™ blister pack.
- Color-coded strip for size identification.
- 5 needles per flat; 100 needle per box.



- Seven package presentations.
- Light-weight aluminum alloy spring handle.
- Innovative and patented package for quick needle access.



For product ID, Prices and specific sizes, please refer to our homepage:
www.houseofacupuncture.dk

For product ID, Prices and specific sizes, please refer to our homepage:
www.houseofacupuncture.dk



NT5: Cluster Pack®



- 5 needles bundled in a regular guide tube.
- Extended medical sealing paper for easy peeling.
- 5 blister packs per flat
- 20 flats; **500 needles per box.**



NT1: Blister Pack



- Individually packaged with SkinSaver™ guide tube.
- 4mm needle length and guide tube differential.
- 5 needles per flat with FingerSaver™ feature.
- 20 flats; 100 needles per box.





- Patented packages are easy to use.
- Lowest package waste per needle in its class.
- Copper handles make the needles easy to locate.
- High density packages with individual needle accessibility.

CT5: Cluster Pack™




- 5 needles bundled in a guide tube.
- 40 flats; 25 needles per flats.
- **1,000 needles per box.**



CTL: Long Needles



- Individually packaged with guide tubes.
- 10 flats; 10 needles per flat.
- 100 needles per box.



For product ID, Prices and specific sizes, please refer to our homepage:
www.houseofacupuncture.dk

For product ID, Prices and specific sizes, please refer to our homepage:
www.houseofacupuncture.dk

7-Star and Intradermal Needles

VINCO® 7-Star Needle Set

- 18 sterile replacement heads, individually packaged and sealed.
- Each head has seven 0.8mm x 3mm needles.
- 6" (15cm) long flexible handle.
- Wrench-style replacement design for maximum user safety.



7-Star needle – Sterile Disposable


- Plastic hammerhead, single sided, black color
- 6" long flexible handle.



Press Tacks

AcuMaster® Press Tacks

- Double-taped for safety
- Individually packaged in blister pack for assured sterility and easy, precise application.
- Made of surgical stainless steel.
- Box of 100.



- APT-01: 0.22mm x 1.5mm
- APT-02: 0.22mm x 1.8mm

VINCO® Press Tack

- **Double-taped** structure for better safety: The needle base is anchored between two layers of surgical adhesive tape.
- Tack will not accidentally fall into ear canals or be left behind when removing it.
- **Extra Comfort:** The base is isolated from the skin and will not cause irritation.
- U.S. - made breathable surgical tape, suitable for sensitive skin.
- 0.22mm dia. x 1.5mm long.
- 10 press tacks on a plastic carrier.



- PT-3415: 100 per box
- 1 - 9 boxes
- 10 - 24 boxes
- 25+ boxes

Ear Pellets

Helio Ear Seeds

- Selected Vaccaria seeds sorted to approximately 2mm size.
- Pre-taped with tan color surgical adhesive tape on plastic carrier. Simply peel and apply!
- Easy-to-use package. Best applied with DP-01 probe (see page 53).
- 100 per box.



- ES-01
1 - 9 boxes
10 - 24 boxes
25+ boxes

Magnetic Ear Pellets


- Approx. 100 gauss with random magnetic pole orientation.
- 1.7mm in diameter with tan colored surgical tape.
- Best applied with DP-01 probe (see page 53).
- Pre-taped.
- 100 pellets per box.



- ES-05
1 - 9 boxes
10 - 24 boxes
25+ boxes

ACU Pellets

- Pre-taped gold or silver plated press pellets.
- 1.6mm in diameter with tan colored round tape.
- 20 pellets per strip; 100 per pack.



- ES-07: Gold
ES-08: Silver

Accu-Patch Press Pellets

- New Improved stainless steel press pellets for the ear and body.
- Similar to Magrain except the 1.2mm diameter stimulating press ball is made from high quality surgical stainless steel on high quality transparent round tape 7mm in diameter.
- Stainless steel.
- 20 pellets per strip. 300 per pack.



- ES-06



Maxim™ - A symbol of quality and value for moxibustion!

MAXIM™ EZY-Lit Smokeless Moxa (Hollow Center)

- This 4.5" x 0.55" Hollow Center EZY-Lit Smokeless moxa sticks are the same composition and size as the familiar Smokeless Moxa Sticks (MX-01).
- No gold paint on the stick.
- Burns at approximately 40 minutes per stick.
- 5 Sticks per box.



- MX-08

MAXIM™ Mini-Smokeless Moxa

- This 0.4" x 3.7" Mini-Smokeless moxa sticks are the same composition as the familiar Smokeless Moxa Sticks (MX-01).
- No Gold paint on the stick.
- Burns at approximately 1.8" per hour.
- Use candle to ignite the roll.
- 10 Sticks per box.



- MX-12

MAXIM™ Smokeless Moxa

- This 0.5" x 4.5" Smokeless moxa sticks are the same composition and size as the familiar Smokeless Moxa Sticks (MX-01).
- No Gold paint on the stick.
- Burns at approximately 1.8" Per hour.
- Use candle to ignite the roll.
- 5 Sticks per box.



- MX-11

Practitioners: To avoid burns, do not leave patients unattended when administering moxa. If the patient complains of discomfort, remove moxa immediately. Moxa is not intended for use on young children or on anyone incapable of detecting the sensation of heat.

MAXIM™ Ultra Pure Jr. Moxa Rolls

- This 4.6" x 0.25" slender tightly packed New Ultra Pure Jr. Moxa Rolls contain a finer grade of moxa.
- Burns at approximately 36 mins per stick.
- 20 Sticks per box.



- MX-07

MAXIM™ Pure Moxa Rolls


- This 0.7" X 8" pure moxa rolls are used for indirect moxa heat treatment.
- The moxa is firmly packed in paper cylinders.
- Pure moxa rolls contain only the leaf of the moxa plant.
- Burns at approximately 4" per hour.
- 10 Sticks per box.



- MX-09

MAXIM™ Stick-On Moxa Cones

- 6mm x 15mm hollow, rigid moxa stick-on cones.
- The base is 0.5mm thick and 10 mm dia.
- The moxa stick is perched on a paper tube to create a thermal buffering zone.
- Burns at approx. 4 mins.
- 180 per box.



- MX-54

Moxa Sticks

Smokeless Moxa

- This 0.5" x 4.5" (1.4cm x 11.4cm) densely packed stick generates intense heat.
- Rigid roll contains moxa and other herbs.
- Very slow burning with a burn rate of approximately 1.8" (4.5cm) per hour.
- Use candle or lighter to ignite the roll.
- 5 Sticks per box.



MX-01

Pure Moxa Rolls

- This 0.6" x 8.0" (1.6cm x 20cm) roll contains loosely packed pure moxa for general moxibustion.
- Burns at approximately 3.5" (9cm) per hour.
- 10 Sticks per box.



MX-02

Nien-Ying Rolls

- 0.6" x 8.0" (1.6cm x 20cm) roll contains tightly packed moxa and herbs.
- Distinguishing aroma produces relaxing effect.
- Burns at approximately 3.4" (8.7cm) per hour.
- 10 Sticks per box.



MX-03

Tai-Yi Medicated Moxa Rolls

- Medicated 0.6" x 8.0" (1.6cm x 20cm) roll contains mugwort (moxa) blended with several natural herbs.
- Use for acute conditions.
- Burns at 4.7" (12cm) per hour.
- 10 Sticks per box.



MX-04

Needle Moxa

Maxim™ Smokeless Needle-Warming Moxa Cones


- Traditional type made especially for Helio, following practitioner's feedback about the excessive weight and complexity of the new wire basket type.
- Lighter, simpler and easier to use than the wire basket type; just place over the needle head and light.
- 450 cones/box.



MX-15

Mini Moxa Rolls

- Used for warming needles. This 0.4"x 0.6" (1.1cm x 1.6cm) firmly packed pure moxa mini rolls yields approximately 6 minutes burning.
- Just fit the needle handle into the hole and squeeze the roll lightly for a tight fit.
- 200 per box.



MX-06

Loose Moxa

Moxa Punk

- MX-20A: Pure Grade
- Yellow green for indirect use.
- 250 gm. per bag



Pure Grade

- MX-20B: Extra Pure
- Fine, yellow green for direct use.
- 50 gm per bag



Extra Pure


- MX-20C: Super Pure
- Golden yellow for direct use.
- 8 gm. per bag



Super Pure

Gold Mountain Ibuki Moxa

- Highest quality moxa available, super pure, extra fine texture, for direct moxibustion.
- Made from carefully selected, superior quality mugwort gathered only in midwinter from Ibuki Mountain when it is at peak quality.
- Pure moxa is less likely to burn the skin because of its high concentration of essential oil.
- Gold Mountain moxa has a volume limit due to its limited gathering and production period.
- Box of 10 gms.



MX-80

Stick-On Moxa Cones

Smokeless Stick-On Moxa Cones

- 6mm x 14mm mini moxa cones
- The base is 0.5mm thick and 10mm dia.
- No smoke
- Pleasant smell
- Mild temperature
- Burns approx. 7 minutes
- 180 per box



MX-53

Tai Ji Mild Stick-On Mini Moxa Cones

- Pure mini moxa cones 4mm x 8.0mm in size.
- Very mild, is ideal for sensitive patients.
- Burning time is approximately four minutes.
- 180 per pack.



MX-56

Skin Shields

- Skin shields are used to reduce heat when using stick-on moxa.
- This skin shields does not prevent burns.
- Do not leave patient unattended during any moxa use.
- 200 per package.



MX-44

Stick-On Moxa Cones

- 7mm x 7mm mini moxa rolls on self adhesive base.
- The base is 4mm thick and 10mm dia with a 3mm center hole.
- Burning time is approximately 3 minutes.
- Use with caution.
- 225 per box.



MX-50

Chosei-Kyu Smokeless Stick On Mini Moxa

- Smokeless and odor free.
- During use no ventilation is required and it does not contaminate fabric with odor.
- Stable heat temperature.
- 100 per box.



MX-79

Sennenkyu Taiyo 24 Self-Heating Moxa

- Self-Heating moxa plasters.
- No fire needed, no smoke created.
- Ideal for patients who cannot tolerate moxa smoke.
- Box of 24 plasters.



MX-78

Moxa Tools

Universal Moxa Extinguisher

- Made of chrome-plated solid brass. 1.75"H x 1.5"dia.
- Extinguishes moxa of various sizes - from 9mm to 16mm in diameter, including smokeless moxa.
- Large heavy base for stability.



MX-36

Moxa Holder/Extinguisher

- Made of chrome-plated solid brass.
- 1.4" high x 1.1" diameter (3.6cm x 2.8cm). Center chamber 0.75" (1.9cm) diameter with a notch on top.
- Can hold a burning moxa roll in the notch or extinguish most moxa rolls in the extinguish chamber in less than 15 seconds.



MX-30

Moxa Extinguisher Dish


- 1.2" diameter moxa roll extinguisher
- Mounted on a wide, stainless steel 4" dia base dish
- Use dish for collecting loose ash from burning moxa



MX-33

Moxage

- Created by merging Japanese craftsmanship with meridian massage and moxibustion.
- Holds heat, reduces smoke, and covers a large area.
- Penetrates deeply to release blocked energy and increase circulation of Qi.
- This system was designed to enhance the healing practices of shiatsu, massage, acupuncture, pediatrics, and veterinary medicine.



MX-34

Needle Caps

- Holds moxa securely on needle handle and ensures no ash can fall onto the skin.
- Air holes allow even burning of moxa.
- Suitable for use only with metal handle needles.
- 10 caps per package.



MX-42

Needle Disposal Containers

Needle Disposal Container

- Two quart, 11"W x 6"H x 4.5"D
- Transparent red
- Used as table top stand alone unit.



NOTE:

The container is environmentally safe. When burned properly at a waste disposal facility, the container emit only carbon dioxide and water. Mailbox-type openings prevent hands from entering the container. A horizontal drop allows for maximum usage of the container volume.

Container accept needles up to 9" (22cm) in overall length

GS-212

Clinical Supplies

Alcohol Prep Pads

- 70% Isopropyl alcohol saturated on a non-woven, bi-folded applicator, packaged in a hermetically sealed foil pouch.
- Sterilized in accordance with U.S.P. sterility testing procedures.
- Used as a topical antiseptic/anti-infective agent prior to administering acupuncture needle insertion.



Medium Size 1.25" x 1.5"

AS-01 box of 200

Large Size 1.75" x 2.0"

AS-02 box of 100

Alcohol Dispenser

- 8 oz. plastic dispensers with white push-down top allows application of alcohol to cotton with one hand.



A. GS-603W

B. GS-603

Diagnostic Tools

Marker Pen

- Ultra fine micro tip, for marking points.
- Gentian violet ink can be removed with alcohol.
- Non-toxic.



MP-01

Pulse Diagnosis Pillow

- A traditional, soft, red velvet pillow as used in China for comfortably supporting the wrist during pulse diagnosis.
- Dimension: 10" x 5" x 3"



GS-818

Cotton Balls

- 500 cotton balls/bag, 100% non-bleached natural cotton.
- Non-sterilized.



GS-601

Medical Supplies

Space Blanket

- Lightweight space blanket keeps your patients feeling warm and safe during treatments.
- Light enough to drape comfortably over inserted needles.
- Size: 56" x 84" (142 x 213 cm)
- Weight: 3 oz. (85g)



GS-805

Cupping Devices

Magnetic Ezy-On Cupping Set


- Magnetic Ezy-On cupping set with 5 plastic cups.
- Including 1 of 1.2", 2 of 1.5", and 2 of 2.3".
- The magnetic probe measures 0.3" (1cm) in diameter and has an intensity of 2500 gauss.
- With vinyl carrying case.
- Made in China.



CP-09

Glass Cupping Jars

- For the classic fire cupping technique.
- Heavy thick glass with smooth mouth.
- Available in three sizes.
- Flat top allows you to rest the jar after cleaning to keep the rim clean.



CP-06B: 1.8" (45mm) in dia
 CP-06C: 2.3" (58mm) in dia
 CP-06D: 2.7" (68mm) in dia

Pistol Hand Pump Cupping Set

- Hand vacuum pump cupping set with 17 plastic cups.
- Including twelve-2.0" (5cm), two-1.75" (4.4cm), one-1.5" (3.8cm), one-1.25" (3.2cm), and one-1.0" (2.5cm) cups in an attractive vinyl case measuring 13.5"x 9.5"x 3.5" (34.5cm x 24cm x 9cm).
- Made in Korea.



CP-03

Replacement cup for CP-03. Specify size.


CP-53

Replacement pump for CP-03

CP-63

Glass Cupping Jars with Rubber Suction Balls

- For the classic cupping technique.
- Heavy thick glass with smooth mouth.
- Available in three sizes.
- Fire unnecessary.
- Made in Germany.



CP-08C: 2.75" (7cm) in dia
 CP-08B: 2" (5cm) in dia
 CP-08A: 1.75" (4.4cm) in dia

You May Purchase from House of Acupuncture with Full Confidence. Manufacturers warranty.

As an ISO certified company, all medical equipment offered by Helio is "Quality Assured" by Helio's rigorous Quality System Management. Helio's warranty policy ensures that your investment is well protected. Our efficient and friendly after-sales service ensures that any need you may have for equipment service is promptly attended. Simply said, Helio is 100% behind all medical equipment offered by Helio.

Helio Equipment Limited Warranty

What does your warranty cover?

- Any defect in material or workmanship.
- Shipping damages.

For how long after purchase?

- Shipping damages reported within five days after receiving the equipment.
- 90 days for labor plus parts.
- One year for parts only.

What will Helio Medical Supplies do?

Full Warranty:

- Repair the defective equipment with no charge for parts and labor within the 90-day period after the purchase date.

Partial Warranty:

- Repair the defective equipment at no charge for parts after 90 days but within the one year period following the date of purchase. Labor for repair will be charged at the prevailing labor rate.
- Helio Medical Supplies, Inc., has the option to provide a new or rebuilt replacement.

How do you get service?

- Call House of Acupuncture at +45 75 81 51 31 to report problems and obtain a return or repair authorization number.
- Customer must prepay all charges when shipping their unit back to House of Acupuncture for repair after the 90-day period from date of purchase. We will refuse the shipment if shipping charges are not prepaid.
- Within the one-year warranty period, we will either repair or replace the unit within 4 weeks after receiving it.

What does your warranty not cover?

- Replacement of owner's manual or operating instructions.
- Installation and set-up service.
- Damages from misuse or neglect.
- When Helio's name plate is removed from the equipment, or when the serial number does not match our record, the warranty is void.
- When the equipment has received unauthorized service or repair.

Helio Repair Service For Firard Series

House of Acupuncture offers a complete line of spare parts (see owner's manual) and repair services for FIRARD Series only. For additional information or pricing call +45 75 81 51 31.

Helio Lamps:

- Qualified technicians will service and repair Firard units. The cost will be \$40 per labor hour plus parts (for a unit purchased from Helio or a Helio's distributor). Typical repair will cost around \$60 plus shipping.
- Turnaround time is about two weeks.
- Please call +45 75 81 51 31 for a repair authorization number or more information prior to shipping your Firard unit.

Non-Helio Lamps:

- Repair for lamps purchased elsewhere: If your original supplier does not provide after-sales support, Helio will charge \$60 per labor hour plus parts. Typical repair will cost around \$90 plus shipping.
- Turnaround time is about two weeks.
- Please call +45 75 81 51 31 for a repair authorization number or more information prior to shipping your unit.

• *Warranty on repaired unit is 30 days from date of invoice.*

Stimulators and Point Detector

Pointer – Excel™ Locator/Stim with Meter

Meter reading increases in relative value from 1 to 500 (at max sensitivity setting) due to a lowering of the electrical resistance of the skin at the acu-point location. Point location is indicated with sound and a pilot light on top of the unit. Point location sensitivity is adjustable.

An effective push button adjustable stimulation offers immediate point treatment. Set "Int" dial to blue zone for micro-amp current output. Includes two probe tips: ear 2 mm (0.08") and body 4 mm (0.16"). Other standard accessories include 9 volt battery, carrying case, and instructions.



EA 901

Self-Calibrating Point Detector

This German-made point detector has an internal memory that memorizes the resistance of the skin. The unit will self-calibrate by the push of a button. Points are located by a light indicator. The unit is 9/16" in diameter x 6 ft" long with a chrome finish. Probe tip is 5/8" long with small tip end for pin-pointing point location. This is extremely useful for ear and hand acupuncture applications. Two 1.5V size-N alkaline batteries are included.



EA 700

POINTER – EXCEL™ SPECIFICATIONS:
 Pulse repetition rate: 1.25 cycles per second pre-set
 Pulse shape: Propriety "M-Wave" (patents applied for)
 Pulse width: 3 Milliseconds pre-set
 Unit Weight: 80 grams
 Unit dimensions: 93 mm x 27 mm (max) x17mm(min)
 Specifications: 500 ohm test load
 Output channels: one
 Frequency: 1 to 16Hz
 Volts & Current: Normal dial: 0 to 22 milliamps
 Blue zone: 0 to 2 microamps
 Polarity reversal: + / - by switch
 Pulse width: 220 microseconds (µS)
 Pulse shape: asymmetric biphasic square wave
 Pulse modes: continuous
 Point detection: meter, sound, and light indication
 Power source: one 9 volt battery
 Size, wt, color: 8.6" x 1.5" x 1", 86gm, beige

Electro-Acupuncture Stimulators

Output Testers
Tests function of channels.

Independently controlled channels
Each channel is adjustable without affecting the output of another.

Safety Feature
E-Strim won't turn on unless all output channels are turned to zero.

Extended Frequency
Adjusts from .5 to 200 Hz

Milliampere Electro-Acupuncture or T.E.N.S compatibility
Selection available

Diagnostic Light

- Check output
- Auto battery check
- Manual battery check
- Clip lead check
- Safety feature

Timer (8C.PRO only)
Machine stays on for the selected amount of time

Calibrated Mixed Mode Frequency
Causes the stimulator to change frequencies every 3 seconds, between two frequencies. Pantheon Research electro-stimulators are unique in their capability to provide adjustment of both frequencies in mixed (dense-disperse) mode.

Back Panel

EA 304, EA 306, EA 308

Electro-Acupuncture Devices

Facial Probe

- Very popular for facial treatments.
- It can also be used on any other part of the body.
- Interchangeable Q-tip/stainless steel ball tips in place of needles.



EA 310

Grounding Probe

Used for those who are afraid of needles, especially children. Similar to the facial probe with interchangeable tips except this has one probe and a hand-held grounding probe.



EA 312

Reusable self-adhesive electrodes

Includes a built-in 2mm wire jack with each electrode.



EA 605: 1,6" x 1,6"
EA 607: 1,5" x 3"

TENS Electrode Lead Wire

3.5mm Jack Plug

Fits EA-304, EA-306 and EA-308 Pantheon & WQ machines, to be used with EA-605 and EA-607 (reusable self-adhesive electrodes).



EA 602

2.5 mm Plug

Fits EA-501, EA-503 and EA-505. This wire also fits EA-2 units allowing you to use electrodes EA-605 and EA-607.



EA 601

2.5 mm Plastic Plug

Fits EA-501



EA 603

Electro-Acupuncture Devices

Model IC-1107+ Electro-Acupuncture Device

Compact, lightweight, sturdy, palm-size, high quality Japanese electro-acupuncture device. A Hi/Low switch on the back of the unit facilitates low voltage setting or Hi voltage setting. This unit is not equipped with a point locator. Accessories include 3 sets of new CE wire leads with rectangular CE jack, carrying case, belt clip, 9 volt battery, and instructions.



EA-905

ITO IC-1107+ SPECIFICATIONS:

Output channels:	3 independent intensity channels
Frequency:	Adjusted with fine frequency dial L – 1 to 20 Hz M – 20 to 150 Hz H – 150 to 500 Hz
Low volts & current:	0 to 10 volts, 0 to 20 milliamps
Size:	3.75" x 2.5" x 1"

SureGrip™ Clips

- An easy push-in plug fit electro-acupuncture equipment
- 24K gold plated for best electric conduction efficiency.
- It weighs only 30mg and is 43% smaller than a miniature alligator clip.
- It is the lightest clip available for electro-acupuncture application.
- The gripping mechanism is a unique innovation.
- There is absolutely no chance for the clip to fall off the needle!



2.5 mm plug (fits EA-2 stimulator)

3.5mm plug (fits most WQ & Pantheon machines)

- EA-201-2: 2 per set
- EA-201-4: 4 per set
- EA-201-8: 8 per set

- WQ-002-2: 2 per set
- WQ-002-4: 4 per set
- WQ-002-8: 8 per set

Electrical Leads

- Electrical leads for electro-acupuncture stimulator.
- 150 cm long leads with two miniature alligator clips on one end and a 3.5mm dia. banana plug on the other end.
- Fits only WQ machines.

WQ-001



TENS

Why Electro-Acupuncture?

Electro-acupuncture (EA) is the stimulation of needles inserted percutaneously with weak alternating electrical currents for therapeutic purposes.

EA greatly enhances acupuncture therapy's effectiveness for a broad range of clinical conditions. Some examples are:

- Accelerates tissue healing of wounds, skin ulcers, etc...
- Facial toning
- Eases acute and chronic pain conditions
- Reduces painful conditions not reacting to manual stimulation
- Decreases substance abuse addictions
- And many others...

TENS, Point Detector and Stimulator

Developed and designed by Helio. This tri-function unit provides the capability of sensitive point detection and electro-acupuncture with accessories that come with the unit. When fitted with TENS accessories, it will provide the function of TENS. EA-2 has a unique three-stage audio indication for point detection: SILENT when the probe is far away from the points, INTERMITTENT sound when the probe is close to a point and CONTINUOUS sound when the probe is on the point. Snap-in battery compartment. Elegant appearance. Soft carrying case for portability. Comprehensive operating manual. A 9V battery is included. Made in China. See page 56 for warranty service.



EA-2

EA-2 TENS SPECIFICATIONS:

Output channels:	Two channels designated as channels A & B
Output Frequency:	Various Modes
Continuous Mode:	1 - 100 Hz, continuously adjustable
Dense Disperse Mode:	Stimulation frequency cycles between 100% and 20% of the selected frequency every eight seconds
Size:	5.25"L x 4"W x 1.5"H

EA-2 Accessories

- EA-200
2.5mm output leads with 2 miniature alligator clips for EA.
- EA-202
2.5mm output leads with banana plugs for point location.
- EA-206
Ground Pole.
- EA-208
Point Locator (Point Probe)



TENS

Micro Current TENS

An advanced micro-current unit that offers both μA and mA stimulation. User-controlled amplitude at fixed frequency. Two selectable wave forms. Recessed dials protect accidental changes of the selected setting. Unit includes two sets of electrode leads and four 1.6" x 1.6" electrodes. Large electrodes may be purchased separately (see page 62). A 9-volt alkaline battery is included. Small size with belt clip for portability.

MICRO CURRENT TENS SPECIFICATIONS:

Channels	Two channels.
Pulse Amplitude	When in μA mode, adjustable 0-900 μA peak into 500W load each channel. When in mA mode, adjustable 0-8mA peak into 500W load each channel.
Pulse Rate	Selectable at 0.6, 10, 30 and 300 Hz
Pulse Width	3, 33, 100 and 166 ms at selected frequency
Wave Form	Asymmetrical bi-phasic square pulse. Selectable square or sloped form
Power Requirement	One 9-volt battery (alkaline or nickel-cadmium rechargeable)
Battery Life	Approximately 100 hours at normal setting
Size and Weight	27 x 62 x 95mm, 130 grams (battery included)



Close up of control panel



EA-505

T-202 TENS

This basic dual-channel TENS unit has full function. User-controlled pulse rate (frequency) and pulse width. Two selectable wave forms in burst and modulation mode. The dial cover protects accidental changes of the selected setting. The unit includes two sets of electrode leads, 4 pieces of reusable rubber electrodes with 12 pieces of disposable adhesive paper (2.75" x 2.5") and one tube of electrode conductivity gel. Accessories can be ordered separately. Self-adhesive electrodes EA-605 and EA-607 can be used with this unit.

T-202 TENS SPECIFICATIONS:

Channels	Two channels.
Pulse Amplitude	0-80 mA each channel, adjustable 500W load.
Pulse Rate	2Hz – 150Hz, adjustable.
Pulse Width	60 μS - 250 μS .
Wave Form	Modified square wave with zero net direct current (DC) components
Operating Mode	Burst and modulation mode.
Power Requirement	One 9-volt battery (alkaline or nickel-cadmium rechargeable).
Battery Life	Approximately 70 hours at normal setting.
Size and Weight	24 x 64 x 95mm, 130 grams (battery included).



EA-507

Close up of control panel



EA-609



EA-742

- EA-507 T-202 TENS set
- EA-609 Disposable adhesive paper for the rubber electrodes
- EA-608 Reusable rubber electrodes, 2 per pack
- EA-742 Electrode conductivity gel, 60ml tube

BioTENS

This basic dual-channel TENS unit offers full function. User-controlled pulse rate (frequency) and pulse width. Three selectable wave forms in continuous, burst and modulation modes. The dial cover protects accidental changes of the selected setting. The unit includes two sets of electrode leads and four 1.6" x 1.6" electrodes. Large electrodes may be purchased separately (See page 62). The 9-volt alkaline battery is included. Small size with belt clip for portability.

BioTENS SPECIFICATIONS:

Channels	Two channels.
Pulse Amplitude	Constant current adjustable 0-80 mA peak into 500W load each channel.
Pulse Rate	2Hz – 120Hz, adjustable.
Pulse Width	40 μS - 260 μS .
Wave Form	Asymmetrical bi-phasic square wave.
Operating Mode	Continuous, burst and pulse width modulation.
Power Requirement	One 9-volt battery (alkaline or nickel-cadmium rechargeable).
Battery Life	Approximately 100 hours at normal setting.
Size and Weight	24 x 64 x 95mm, 130 grams (battery included).



Close up of control panel



EA-501



The World's #1 Brand in Massage

AVILA II™ Massage Chair



Avila II™ Massage Chair Complete Package:
Includes carry case with wheels, sternum pad and instructional video.

- Weight: 22 lbs
- Standing: 34"L x 18 1/2"W x 47"H
- Folded: 44"L x 20"W x 12 3/4"H

ET-105




AVALON™ Massage Table



The new Avalon™ table combines value and unparalleled performance for the massage professional or serious student. Feature after feature of tables costing much more - are standard with our Avalon. It's a workhorse table that's handsome and durable, too - with a lifetime limited warranty to give you peace of mind. Avalon table is the ultimate combination of honest value, intelligent design and true craftsmanship.

- Weight range: 32-33 lbs
- Width: 30"
- Height range: 23" to 33"
- Cradle Lock Cable System
- Lifetime Limited Warranty
- Premium maple hardwood and Russian birch decking



Avalon Table Only:

ET-053: 32 lbs, 30"x73"

Avalon™ Student Package:

Includes table, deluxe adjustable headrest, and single-pocket carry case.

ET-053P: 32 lbs, 30"x73"




The World's #1 Brand in Massage

LUNA™ Massage Table



Superlite cutting edge frame design. Full-sized, lightweight, and comfortable.

- Reiki ready lightweight aluminum frame, patent pending
- Weight range: 26-29 lbs
- Width: 30"
- Height range: 24" to 34"
- Cradle lock cable system
- Lifetime Limited Warranty



Luna™ Table Only

ET-178: 26 lbs, 30"x73"

Luna™ Silver Package

Includes table, and NEW! Caress adjustable face cradle and standard single-pocket carry case.

ET-178P: 26 lbs, 30"x73"

HARMONY 3™ Massage Table



The Harmony 3™ is our most economically priced full size massage table. Built with Russian birch decking and maple frame. The value package can't be beat.

- Weight: 29.5 lbs
- Width: 29"
- Height range: 23" to 34"
- Cradle lock cable system
- Lifetime Limited Warranty



Harmony 3™ Table Only

ET-060: 31 lbs, 29"x73"

Harmony 3™ Complete Package

Includes table, deluxe headrest-black, and black carry case.

ET-675: 31 lbs, 29"x73"





Massage Table Accessories

Deluxe Adjustable Headrest

ET-304: Black Frame
ET-303: Cream Frame



Bolsters: Jumbo Round

6"x26"
ET-357-A: Agate
ET-357-T: Teal
ET-357-B: Black
ET-357-P: Purple



Bolsters: Peak

7"x8 1/2"x24"
ET-360-A: Agate
ET-360-T: Teal
ET-360-B: Black
ET-360-P: Purple



Bolsters: Half Jumbo

4.5"x9"x29"
ET-355-A: Agate
ET-355-T: Teal
ET-355-B: Black
ET-355-P: Purple



Bolsters: Fluffy

8"x26"
ET-352-A: Agate
ET-352-T: Teal
ET-352-B: Black
ET-352-P: Purple



Bolsters: Full Round

6"x29"
ET-353-A: Agate
ET-353-T: Teal
ET-353-B: Black
ET-353-P: Purple



Human Body

Acupuncture Meridian Point Models

HM-01, HM-02, and HM-03 models are made of soft white vinyl, while HM-07 is made of a solid fiberglass-type material. All have 360 points and 14 meridian in both Chinese and standard English designations. Also shown are 80 extra non-meridian points. Wooden display stand included.



HM-01: 18" (46cm) male model
HM-02: 18" (46cm) female model
HM-03: 22" (56cm) male model
HM-07: 27" (69cm) male model

Shipping: 12 lbs for HM-07. Shipping charges at actual shipping weight.

Human Head Model


Soft vinyl 8.5"(22cm) high head model on wooden base. Meridian points are shown on head, face and scalp. All points are identified by both Chinese and standard English designations.



HM-05

Acupuncture Book Charts

Set of 4 acupuncture charts. Each measures 13"x 26" (33cm x 66cm) in a bound 9.5" x 13" (24cm x 33cm) book. Very high quality full-color printing on heavy bond paper. Alpha-numeric and Chinese names are both presented. Printed in Hong Kong.



BC-102

Illustrated Acupuncture Chart

A 24"x 40" (71cm x 114cm) chart by Michael Spatuzzi. Shown are the phase points, accumulation, hoary, connecting and source, tonification and sedation points of the 12 main meridians as well as the window of the sky, influential points and points of the four seas. Printed in the USA.




BC-109

Human Body

Acupuncture Wall Chart Set


Set of 4 acupuncture charts. Each measures 21"x 44" (53cm x 112cm). Very high quality full-color printing on heavy bond paper. Alpha-numeric and Chinese names are both presented. Extremely high quality prints. Printed in Hong Kong.

BC-101



Trigger Point Charts

- Clearly displays locations and pain patterns of trigger points in all major muscles.
- Each muscle shown on a separate panel.
- Vital information for anyone who works with pain.
- Designed for the professional.
- A collaboration between a clinician with 20 years experience and a graphic artist specializing in anatomical illustration.



Set of I and II
BC-120

Ear

Ear Models

HS-61: 9" / 23cm

Made of peach-colored soft plastic with points numbered from 1 to 41. A wooden stand is included.

HS-62: 5" / 13cm

Made of peach-colored soft plastic with points numbered from 1 to 41

HS-64: 5" / 13cm

Soft vinyl with light skin tone 5" long ear acupuncture model indicating organ names.

HS-65: 7"(17cm)

Made of peach-colored soft plastic with points indicating organ names in English.

HM-63A: Right Ear

HM-63B: Left Ear


Made of soft silicone rubber for realistic skin feel. Made in Germany.



Auricular Therapy Reflex Desk Reference Cards

A two-sided 8.5" x 11" Quick Reference Card for keeping on the counter in your treatment room. Contains the same information as the large Wall Chart, with Phases 1 and 4 on the front and Phases 2 and 3 on the reverse. Card is constructed of heavy glossy card stock for ease of use.


BC-113



Auricular Therapy Reflex Chart

Photo-finish quality 24"x 36" developed by Bryan L. Frank, M.D. and Nader E. Soliman, M.D. Suitable for wall display.

BC-111



Ear

Ear Reflex Chart

By Terry Oleson

This photo-finish quality 24"x 38" (61cm x 97cm) color chart shows surface and hidden views of the ear for more than 200 ear points in both Chinese and French auricular medicine. Points are designated by numbers in English and Chinese. Printed in the USA.


BC-105



Auriculotherapy Manual, 3rd Edition

By Terry Oleson, Ph.D.
BC-500

This handbook combines the best work of Chinese and French systems of ear acupuncture-auriculotherapy. The ear reflex points presented in this handbook are based upon extensive research conducted at the UCLA Pain Management Center. This resource provides a standardized auricular nomenclature for designating different anatomical zones of the ear. 260 pages




Pocketbook Atlas of Auricular Therapy and Auricular Medicine



By Bryan L. Frank, M.D. and Nader E. Soliman, M.D.

Same illustrations as the highly acclaimed Atlas of Auricular Therapy and Auricular Medicine in a new pocket-sized (4.25"x 5.5") Atlas! Now you can have the leading advanced Auricular Therapy Atlas and the convenience and portability of the pocket-sized version.


BC-116



Atlas of Auricular Therapy & Auricular Medicine

The Atlas provides a perfect complement to the Reference Charts, and is designed on a map design rather than an anatomic design, to allow identification of structures that cannot be drawn on the anatomic charts. Various reference zones from the works of Nogier are included, as well as treatment indications from the authors' own experience. The Atlas is constructed with a spiral spine for ease of use and is on 24-pound heavy paper. It contains 150 illustrations, 138 of which are in full color.

BC-115



Hand and Foot

Hand Models


HS-70
Small, five inches

HS-72
Large, eight inches with a wooden stand



Hand Therapy Charts

Beautiful anatomic rendering with 2-sided, full color illustrations. These 11"x17" charts use conventional acupuncture channel numbering rather than confusing numbering systems often seen. Command Points and Influential Point reference are given.



BC-112

Foot Model

Made of soft vinyl 5.5" foot acupuncture model showing acupuncture points around the ankle, top and bottom of the foot.

HS-75



More Than One Way To Serve You!

House of Acupuncture is the sole importer of Helio® products for Denmark.

Many fine products such as AcuGlide®, AcuMaster®, VINCO® MicroClean®,

Viva™, and C&G™ acupuncture needles, FIRARD II™ and FIRARD-3™ far

infrared therapeutic devices, EA-2 and Pantheon multiple-function electro-

acupuncture/TENS stimulator and Maxim™ Moxibustion products to name a few.

Danish importer of Helio® products:

House Of Acupuncture Aps

Postbox 1029

7120 Vejle Ø

Denmark

Phone. +45 75 81 51 31

Fax. +45 75 72 56 16

E-mail. shop@houseofacupuncture.dk

Website. www.houseofacupuncture.dk

House of Acupuncture

danish distributor of



Have you visited our NEW site?

Check out our new look and these new features:

- Order online and track order status
- Account maintenance, including password, email and contact info change
- "What's New" section highlights new products, seasonal promotion and upcoming events
- "Customer Support" if you have any questions or to request free catalog and samples

www.houseofacupuncture.dk

